# Lab 1 Motion diagrams

# Rectangle Repulsed Researchers

Julian Barossi, Liam Gilligan, Stephanie L'Heureux

September 6, 2023

Part 3: Matching Position vs. Time and Velocity vs. Time Graphs

## Graph 1




Figure 1: Graph 1

## Graph discription

Graph one depicts a the position vs time graph Graph 1 shows a graph of a collision. The object's distance from the sensor is increasing. The slope of the position is velocity so we can also say that the velocity is positive and constant. Then at  $t \approx 2.4s$  its position begins decreasing at a constant, negative velocity.

$$\frac{d}{dx}[s(t)] = v_{x_1} = 0.35m/sV_{x_2} = -0.35m/s \tag{1}$$

Setup

Our setup consisted as a flat track with a magnetic bumper 60cm from the the *PASCO Universal 850 interface and Motion Sensor II*. We pushed and released the cart towards the bumper. We consider friction and drag negligible.

#### Graph 2:

Graph discription:

Graph 2 depicts an object speeding up at a constant rate. We know this because the object's position is increasing and the velocity is positive and increasing.

Setup:

Our setup consisted of a track at an downward incline. We released the cart from the top of the track and let it slide to the bottom. We consider friction and drag negligible.

#### Graph 3:

Graph discription:

Graph 3 depicts an object in freefall that hits the ground and stays there. We know this because the object's position starts from an increased height

Setup:

Our setup consisted of a the sensor positioned updward. We dropped a large beach ball.

### Graph 4:

Graph discription:

Graph 4 depicts an object moving down an incline. It comes in contact with something at the bottom and changes direction. As it changes direction it slows down. Then at t=2.6s its velocity becomes negative but the acceleration staus positive

where it is slowing down

Setup:

Our setup consisted of a the sensor positioned at the biotton of an incline plane. We released the cart from the top

Graph 5:

Graph discription:

Graph 5 depicts an object that starts from rest then is that ts pushed then slows to a stop. We know this because the graph depicts the initial velocity as 0m/s, then the velocity rapidluy increases and the acceleration is positive. It the acceleration goes to zero and he velocity is constant

Setup:

Our setup consisted a flat plane and w pushed the cart away from the sensor

Graph 6:

Graph discription:

Graph 6 depicts an object which is pushed and released up a hill. Eventually gravity slows it down and it returns down the hill

Setup:

Our setup consisted of a incline plane. We gave the cart a push and let it go up the hill then return down

Graph 6:

Graph discription:

Graph 6 depicts an object moving up and swown on a spring. Setup:

Our setup consisted of an downward incline. The bumper was positioned 60cm away from the sensor.